


ROSA DONQ2 Quantitative Test for Feed & Grain

2
minutes

Product Overview



The **ROSA® DONQ2 Quantitative Test (DONQ2)** delivers rapid and accurate results using ROSA (Rapid One Step Assay) lateral flow technology. Quantitative results are available in a simple, two-minute assay using water as the extraction solvent to detect DON (vomitoxin) in feed and grain.

Benefits of the DONQ2 Test include:

- ◆ Range: 0 to 5.4 ppm, 2 to 25 ppm, 5 to 60 ppm
- ◆ Results in 2 minutes
- ◆ GIPSA Approval
- ◆ Negative samples may be disposed as normal waste
- ◆ Awarded 5 year GIPSA contract to monitor DON (vomitoxin) in grains

Run on the Charm EZ-M System:

- ◆ High volume output with batching using multiple quad incubators
- ◆ Standard curve built into reader reducing testing cost by eliminating end-user calibration
- ◆ Printable or downloadable results to any computer or LIMS system

Simple Procedure

Add Sample.

Incubate.

Read results in less than 5 seconds!

For batch processing, incubate multiple strips simultaneously, then read results in the Charm EZ-M system.

GIPSA^A Approved Commodities

Barley	Brewer's Rice	Brown Rice	Buckwheat	Corn
Corn Bran	Corn Germ Meal	Corn Gluten Feed	Corn Gluten Meal	DDGS
Hominy	Malted Barley	Milled Rice	Oats	Rapeseed Meal
Rice Bran	Rough Rice	Rye	Sorghum	Soybean Meal
Triticale	Wheat	Wheat Bran	Wheat Flour	Wheat Midds
Wheat Red Dog	Wheat Screenings			

^AUSDA - Grain Inspection, Packers and Stockyard Administration - Federal Grain Inspection Services

Charm Validated Commodities

Palm Kernel Meal				
------------------	--	--	--	--

Ordering Info

Order Codes

LF-DONQ2-20K
LF-DONQ2-100K
LF-DONQ2-500K

Each Kit Includes

Test Strips, Positive Control, Dilution Buffer

See Operator's Manual for additional kits and order codes



659 Andover Street, Lawrence, MA 01843-1032
T +1.978.687.9200 | F +1.978.687.9216
E info@charm.com | www.charm.com

© 2017 Charm Sciences, Inc. Charm, Charm EZ, and ROSA are registered trademarks, and test strip product packaging is trade dress, of Charm Sciences, Inc.

MRK-933-005 Apr-2017